

**SYRACUSE MEDICAL DEVICES, INC.**

214 Hurlburt Road, Syracuse, NY 13224-1821

Tel.. 315-449-0657 FAX 315-449-0756

Email: [syrmed@twcny.rr.com](mailto:syrmed@twcny.rr.com)

June 26, 2024

The vaginal dilators that Syracuse Medical Devices manufactures are a hard plastic material. Our dilators are composed of 100% high-density polyethylene and have zero latex content. For sterilization you may use a gas autoclave (ETO) and a temperature range not to exceed 200 degrees Fahrenheit.

Syracuse Medical Devices Vaginal Dilators are considered disposable as marketed, single patient use, and are priced accordingly.

Patient instructions for cleaning are to use hot water and a mild detergent.